

Introduction

RNA capping in humans is essential for mRNA splicing, protection of RNA from 5'-exonucleases in the cytoplasm and targeting of mRNAs to the ribosome. Human RNA guanine-7 methyltransferase (RNMT), Cap Methyltransferase 1 (CMTR1), and Cap Methyltransferase 2 (CMTR2) are RNA methyltransferases involved in RNA capping, play significant roles in the proliferation and differentiation of embryonic stem cells, and have also been implicated in cancer.

A study by Taherian et al. provided a comprehensive understanding of the substrate specificities of human RNA capping methyltransferases, essential for the development of future anticancer therapeutics and assessment of antiviral therapeutics' selectivity. The ability to synthesize highly modified and structurally diverse RNA molecules is crucial for advancing gene editing technologies and nucleic acid-based therapeutics.

Here, we utilize an integrated strategy that combines solid-phase synthesis with enzymatic methods to achieve high-yield production of functional RNA constructs, including 5'-capped RNA (TMG- and m⁷G).

In eukaryotes, capping is executed by a series of highly coordinated enzymatic steps, typically recruited by the phosphorylated C-terminal domain (CTD) of RNA Polymerase II:

1. RNA Triphosphatase:

Removes the γ -phosphate from the 5'-triphosphate end of the nascent pre-mRNA, leaving a 5'-diphosphate end (5').

2. RNA Guanylyltransferase:

Transfers a GMP moiety from a GTP molecule to the 5'-diphosphate, forming the unique 5'5' – triphosphate bridge. This results in an unmethylated G(5')ppp(5')N structure.

3. RNA (Guanine-N7)-Methyltransferase:

Methylates the guanine base at the N7 position, using S-adenosylmethionine (SAM) as the methyl donor to yield the definitive m⁷G cap.


Depending on subsequent methylations of the adjacent transcribed nucleotides, caps are classified into different variations:

Cap 0: Only the terminal guanine is methylated (m⁷GpppN). This is typical in lower eukaryotes like yeast.

Cap 1: The 2'-hydroxy (2) group of the first transcribed ribose nucleotide is additionally methylated by a 2'-O-methyltransferase (m⁷GpppN_m). This is the predominant form in higher eukaryotes and serves as a crucial marker to prevent the innate immune system from recognizing the mRNA as foreign.

Cap 2: The 2'-OH groups of both the first and second transcribed nucleotides are methylated (m⁷GpppN_mN_m).

RNA Caps





Synthesis Strategy

- [1] Solid phase synthesis of PRE-oligonucleotides
- [2] Purification by HPLC
- [3] Quality Control
- [4] Addition of functional groups using enzymes
- [5] Purification by HPLC
- [6] Final Quality Control

Synthesis of modified RNA with a 2,2,7-trimethylguanosine 5'-cap structures (TMG-cap)

The 2,2,7-trimethylguanosine (TMG) 5' cap is an important modification that enhances RNA stability, facilitates nuclear transport, and promotes translation, particularly in snRNAs and certain viral RNAs.

Name	Sequence (5' – 3')	Length	MW
RNA_oligo2	[TMGppp]mAmGmUmUmAGCTGGATTACGGGGTTCATTAGTCCAGTACITTCATAG	42	13,644.68
U7 snRNA	[TMGppp]Cra rGUg rUUA rCrAg rUUC rUUU rUAAG rAAU rUUg rUCU rAGrU rAGG rUUrU rUCrU rGGrC rUUU rUUUA rCrCG rGAaA rAGrC rCrC rU	63	20,544.96




ESI-MS

HPLC

The 2,2,7-trimethylguanosine (TMG) cap and capped synthetic U7 snRNA.

Synthesis of modified RNA with a Guanosine 5'-cap structure (Gppp-cap)

Sequence (5' – 3')	Calculate MW
[Gppp][mA]rCr rCrCrC rCrCrC rCrC[3-Biotin]	3970.51




ESI-MS Result (A) and HPLC profile (B) of RNA_oligo1 synthesized using a combined chemical and enzymatic approach.

Synthesis of a highly modified RNA with 7-methylguanosine 5' cap structures (m⁷G-cap)

Large-scale synthesis of m⁷G-capped RNA with high purity for the production of high-quality capped RNAs for diverse downstream applications.

Name	Sequence (5' – 3')	MW
RNA_oligo2	[M7Gppp]rGrArCrUrUrCrGrArCrGrCrUrCrArUrCrGrArGrUrArUrArGrA[AmC6]-[5-6ROX]	9,197.73



ESI-MS analysis (A) and HPLC profile (B) of m7GpppG-RNA oligonucleotides synthesized using a combined chemical and enzymatic approach.

Conclusion

We show that an integrated synthesis strategy combining solid-phase synthesis with enzymatic modifications enables the production of high-quality, structurally diverse nucleic acids. Our optimized workflow enables the efficient generation of highly modified RNA molecules, including 5'-capped (TMG- and m⁷G-) RNAs, ensure enhanced stability and functionality.

By integrating chemical synthesis precision with enzymatic efficiency, our approach overcomes the limitations of traditional methods such as in vitro transcription. This scalable, reproducible, and high-yield synthesis platform facilitates advancements in gene-editing technologies, RNA therapeutics, and molecular diagnostics.

References

- Chan SH, Molé CN, Nye D, et al. Biochemical characterization of mRNA capping enzyme from Faustovirus. *RNA*. 2023;29(11):1803–1817. doi:10.1261/ma.079738.123.
- Ignatovich AV, Iguchi JA, Kore AR, Ho CK. Trypanosome mRNA recapping is triggered by hypermethylation originating from cap 4. *Nucleic Acids Res*. 2024;52(17):10645–10653. doi:10.1093/nar/gkac614.
- Molé CN, Chan SH, et al. Kinetics study reveals interconnected relationships of component activities in trifunctional RNA capping enzymes. *J Biol Chem*. 2026;302(8):113237.
- Podvalnaya N, Bronkhorst AW, Lichtenberger R, et al. piRNA processing by a trimeric Schlafen-domain nuclease. *Nature*. 2023;622(7982):402–409. doi:10.1038/s41586-023-06588-2.
- Taherian F, Mabanglo MF, Hajian T, et al. Defining substrate specificities of human RNA capping methyltransferases through quantitative assessment of independent yet cooperative activities. *Protein Sci*. 2026;35(1):e70428. doi:10.1002/pro.70428.
- Wang J, Chew BLA, Lai Y, et al. Quantifying the RNA cap epitranscriptome reveals novel caps in cellular and viral RNA. *Nucleic Acids Res*. 2019;47(20):e130. doi:10.1093/nar/gkz751
- Wolf EJ, Grünberg S, Dai N, et al. Human RNase 4 improves mRNA sequence characterization by LC-MS/MS. *Nucleic Acids Res*. 2022;50(18):e106. doi:10.1093/nar/gkac632.