

# Sensitive Detection of BRAF V600E Mutation by BNA Clamping Real-time PCR

Xiaoyun Liu<sup>1</sup>, Houquan Dai<sup>1</sup>, Aaron Castro<sup>1</sup>, Leticia Loredo<sup>1</sup>, Yuewei Zhao<sup>1</sup>, Sung-Kun Kim<sup>2</sup>, Austin Dinkel<sup>2</sup>, Miguel Castro<sup>1</sup>

<sup>1</sup>Biosynthesis Inc, 612 E. Main St, Lewisville, TX 75057, USA

<sup>2</sup>Department of Natural Sciences, Northeastern State University, 611 N. Grand Ave, Tahlequah, OK 74464, USA



## Purpose

We report a sensitive method for the rapid detection of the 1799T>A (V600E) conversion on BRAF gene, using 2'-O,4'-aminoethylene bridged nucleic acid (BNA or BNA<sup>NC</sup>). We have designed, synthesized, and investigated several BNA clamps for their superior binding ability to perfectly matched DNA templates while discriminating the BRAF mutant from the wild-type gene. We have also produced a fluorescence-labelled BNA probe for specific detection of the BRAF-V600E mutation. We have found that the BNA clamp-probe combination is able to detect mutants at abundance levels as low as 0.1%, indicating an improved sensitivity for the diagnosis of gene mutations.





Figure 1. The structures of DNA and synthetic analogues.



Figure 2. Illustration of the B-raf protein kinase and the location of the V600E mutation.

## Methods

- BNA clamp and probe-based real-time PCR



## Key Findings


- Incorporation of BNA bases significantly increases melting temperatures of the oligonucleotides.
- The inhibitory effect ( $\Delta Ct$ ) of the BNA clamps on the amplification of the WT gene, directly correlates to their  $T_m$  values, while BNA clamp-1412R shows the highest efficacy ( $\Delta Ct > 12$ ).
- The BNA probe is able to specifically detect the BRAF-600E mutant. It is able to detect levels as low as 10 cps of the mutant gene.
- In combination with BNA Clamp-1412R, the BNA probe is able to detect a mutation level as low as 0.1% in a background of the wild-type gene.**
- BNA clamp-1412R and BNA probe is able to specifically detect the BRAF-V600E mutation present in genomic DNA.

## Results

Table 1. BNA-clamping oligonucleotides and  $T_m$  values for synthetic ssDNA 14-mer templates of the BRAF 600V WT and 600E mutant genes

Name	Length	No. of BNA	$T_m$		
			WT	600E	$\Delta T_m$
1408F	14	8	68.5	54.5	14.0
1411F	14	11	80.1	68.2	11.9
1412F	14	12	80.7	65.1	15.6
DNA-F control	14	0	48.5	37.5	11.0
<b>1412R</b>	14	12	83.6	66	17.4
DNA-R control	14	0	46.5	37.5	9.0

Figure 3. Real-time PCR amplification of BRAF genes in the absence or presence of the BNA clamps. A total  $10^3$  copies of plasmids were used in each assay.



## Results

Figure 4. Direct correlations of the melting temperatures of BNA clamps with the inhibitory effect on the PCR amplification of the WT gene




Figure 5. Real-time PCR detection of the BRAF 600E mutant by the BNA and universal probes





Figure 6. Ultra sensitive detection of the BRAF-V600E mutant by BNA clamp-1412R and probe.



## Results

Figure 7. Real-time PCR detection of the BRAF-V600E mutation in genomic DNA (gDNA) samples by BRAF BNA clamp-1412R and BNA probe.




Table 2. Real-time PCR detection by the BNA clamp-1412R and BNA probe, of the BRAF-V600E mutant present in the genomic DNA samples

Samples	BRAF Status	Assay	Ct
Human gDNA	600V (WT)	Total	32.4 ± 0.3
		Mutant	> 50
SW480	600V (WT)	Total	26.7 ± 0.1
		Mutant	> 50
HT29	600E (heterozygous)	Total	26.2 ± 0.0
		Mutant	30.4 ± 0.1
A375	600E (homozygous)	Total	28.2 ± 0.2
		Mutant	30.6 ± 0.1
NTC	N/A	Total	>50
		Mutant	>50

## Future Directions

The BNA clamp-1412R and probe assay will be tested on DNA extracted from formalin-fixed, paraffin-embedded (FFPE) human tissues and freshly isolated tumor samples.

## References

- Hirama T, Shiono A, Egashira H, et al. PCR-Based Rapid Identification System Using Bridged Nucleic Acids for Detection of Clarithromycin-Resistant Mycobacterium avium-M. intracellulare Complex Isolates. J Clin Microbiol. 2016;54(3):699-704.
- Morishita S, Takahashi K, Araki M, et al. Melting Curve Analysis after T Allele Enrichment (MelcaTle) as a Highly Sensitive and Reliable Method for Detecting the JAK2V617F Mutation. PLoS One. 2015;10(3):e0122003.