

Truly Total RNA Purification Kit

Excellent sensitivity, Superior microRNA recovery, Extraction from a single cell

FEATURES AND BENEFITS

- ▶ No phenol chloroform extractions
- ▶ No need for carrier RNA
- ▶ 3 Easy Steps: Bind, Wash, Elute
- ▶ Isolate total RNA from a broad input source from animal cells, tissue, bacteria, yeast, fungi, plants, viruses and various bodily fluid
- ▶ Isolate a diversity of RNA species from mRNA to microRNA & siRNA
- ▶ Fast and easy processing based on spin-column format in 20 minutes
- ▶ Isolate total RNA from very small samples

Purification is based on spin column chromatography using BSI proprietary resin as the separation matrix.

We have a kit for that...

- ▶ DNA Purification Kits
- ▶ Protein Purification Kits
- ▶ Urine Sample Preparation Kits
- ▶ Blood/Plasma/Serum Sample Preparation Kits
- ▶ Molecular Diagnostics
- ▶ PCR reagents



BIO SYN.COM

Service across Genomic, Proteomic, and Cell Biology

BIO SYNTHESIS
COMMITTED TO BIOMIC RESEARCH

Toll Free: 800.227.0627 | 972.420.8505

info@biosyn.com

612 E. Main Street, Lewisville, TX 75057 USA